

## HES1 antibody

### Product Information


Catalog No.:	FNab10698
Conjugation:	None
Size:	100µg
Form:	liquid
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	HES1  BHLHB39  HL  HRY
Uniprot ID :	Q14469

### Application

Reactivity:	Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:1000



Mouse eye lysates were subjected to SDS PAGE followed by western blot with FNab10698(HES1 antibody) at dilution of 1:500