

---

## Goat anti- monkey IgM

### Product Information

Catalog No.:	FNab10700
Size:	1mg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Goat
Clonality:	polyclonal
IsoType:	IgG
Storage:	10mMPBS,pH 7.3, 4°C for 12 months.

### Immunogen information

Immunogen:	monkey IgM mu-chains
------------	----------------------

### Application

Reactivity:	Monkey
Tested Application:	ELISA, WB

**Note: This product is for research use only, not for diagnostic or therapeutic use.**