

## P4HA1 antibody

### Product Information

Catalog No.:	FNab10712
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Catalyzes the post-translational formation of 4-hydroxyproline in-Xaa-Pro-Gly-sequences in collagens and other proteins.

### Immunogen information


Immunogen:	prolyl 4-hydroxylase, alpha polypeptide I
Synonyms:	Prolyl 4-hydroxylase subunit alpha-1 (4-PH alpha-1) Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-1 P4HA1 P4HA
Observed MW:	61 kDa
Uniprot ID :	P13674

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10712(P4HA1 antibody) at dilution of 1:100



A431 cells were subjected to SDS PAGE followed by western blot with FNab10712(P4HA1 antibody) at dilution of 1:600