

## DAB2 antibody

### Product Information


Catalog No.:	FNab10718
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	disabled homolog 2, mitogen-responsive phosphoprotein(Drosophila)
Synonyms:	Disabled homolog 2 Adaptor molecule disabled-2 Differentially expressed in ovarian carcinoma 2 (DOC-2) Differentially-expressed protein 2 DAB2 DOC2
Observed MW:	58 kDa, 96 kDa
Uniprot ID :	P98082

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate tissue slide using FNab10718(DAB2 antibody) at dilution of 1:100



HeLa cell lysates were subjected to SDS PAGE followed by western blot with FNab10718(DAB2 antibody) at dilution of 1:1000