

ZNF140 antibody

Product Information

Catalog No.:	FNab10746
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation as a repressor.

Immunogen information


Immunogen:	zinc finger protein 140
Synonyms:	Zinc finger protein 140 ZNF140
Observed MW:	41 kDa, 53 kDa
Uniprot ID :	P52738

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab10746(ZNF140 antibody) at dilution of 1:100



U87 cell lysates were subjected to SDS PAGE followed by western blot with FNab10746(ZNF140 antibody) at dilution of 1:1000