

CHGA antibody

Product Information


Catalog No.:	FNab10750
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	chromogranin A
Synonyms:	Chromogranin-A (CgA) Pituitary secretory protein I (SP-I) Vasostatin-1 Alternative names: Vasostatin I Vasostatin-2 Alternative names: Vasostatin II EA-92 ES-43 Pancreastatin SS-18 WA-8 WE-14 LF- 19 Catestatin Alternative names: SL21 AL-11 GV-19 GR-44 ER-37 GE- 25 Serpinin-RRG Serpinin p-Glu serpinin precursor CHGA
Observed MW:	50-100 kDa
Uniprot ID :	P10645

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10750(CHGA antibody) at dilution of 1:100



LO2 cell lysates were subjected to SDS PAGE followed by western blot with FNab10750(CHGA antibody) at dilution of 1:1000