

GOLM1 antibody

Product Information


Catalog No.:	FNab10755
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	golgi membrane protein 1
Synonyms:	Golgi membrane protein 1 Golgi membrane protein GP73 Golgi phosphoprotein 2 GOLM1 C9orf155 GOLPH2
Observed MW:	73 kDa
Uniprot ID :	Q8NBJ4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cell lysates were subjected to SDS PAGE followed by western blot with FNab10755(GP73 antibody) at dilution of 1:1000