

FLCN antibody

Product Information

Catalog No.:	FNab10758
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be a tumor suppressor. May be involved in energy and/or nutrient sensing through the AMPK and mTOR signaling pathways. May regulate phosphorylation of RPS6KB1.

Immunogen information

Immunogen:	folliculin
Synonyms:	Folliculin BHD skin lesion fibrofolliculoma protein Birt-Hogg-Dube syndrome protein FLCN BHD
Observed MW:	75 kDa
Uniprot ID :	Q8NFG4

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A431 cell lysates were subjected to SDS PAGE followed by western blot with FNab10758(FLCN antibody) at dilution of 1:1000