

CEP350 antibody

Product Information


Catalog No.:	FNab10770
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	centrosomal protein 350kDa
Synonyms:	Centrosome-associated protein 350 (Cep350) Centrosome-associated protein of 350 kDa CEP350 CAP350 KIAA0480
Observed MW:	200 kDa
Uniprot ID :	Q5VT06

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded rat brain slide using FNab10770(CEP350 Antibody) at dilution of 1:50