

PRF1 antibody

Product Information


Catalog No.:	FNab10778
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	perforin 1(pore forming protein)
Synonyms:	Perforin-1 (P1) Cytolysin Lymphocyte pore-forming protein (PFP) PRF1 PFP
Observed MW:	None
Uniprot ID :	P14222

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat lung tissue slide using FNab10778(perforin Antibody) at dilution of 1:100