

ITFG1 antibody

Product Information

Catalog No.:	FNab10782
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Modulator of T-cell function. Has a protective effect in graft versus host disease model(By similarity).

Immunogen information

Immunogen:	integrin alpha FG-GAP repeat containing 1
Synonyms:	T-cell immunomodulatory protein (Protein TIP) Integrin-alpha FG-GAP repeat-containing protein 1 Linkin ITFG1 LNKN-1 TIP
Observed MW:	None
Uniprot ID :	Q8TB96

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10782(ITFG1 Antibody) at dilution of 1:100