

## GCG antibody

### Product Information


Catalog No.:	FNab10785
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	glucagon
Synonyms:	Pro-glucagon Glicentin Glicentin-related polypeptide (GRPP) Oxyntomodulin (OXM, OXY) Glucagon Glucagon-like peptide 1 (GLP-1) Alternative names: Incretin hormone Glucagon-like peptide 1(7-37) (GLP-1(7-37)) Glucagon-like peptide 1(7-36) (GLP-1(7-36)) Glucagon-like peptide 2 (GLP-2) GCG
Observed MW:	None
Uniprot ID :	P01275

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded human pancreas using FNaB10785(GCG antibody) at dilution of 1:200