

CD36 antibody

Product Information


Catalog No.:	FNab10787
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	CD36 molecule(thrombospondin receptor)
Synonyms:	Platelet glycoprotein 4 Fatty acid translocase (FAT) Glycoprotein IIIb (GPIIIB) Leukocyte differentiation antigen CD36 PAS IV PAS-4 Platelet collagen receptor Platelet glycoprotein IV (GPIV) Thrombospondin receptor CD36 GP3B GP4
Observed MW:	32-53 kDa
Uniprot ID :	P16671

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNaB10787(CD36 antibody) at dilution of 1:100



HepG2 cell lysates were subjected to SDS PAGE followed by western blot with FNaB10787(CD36 antibody) at dilution of 1:1000