

SPINT1 antibody

Product Information

Catalog No.:	FNab10792
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


The protein encoded by this gene is a member of the Kunitz family of serine protease inhibitors. The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms.

Immunogen information


Immunogen:	serine peptidase inhibitor, Kunitz type 1
Synonyms:	Kunitz-type protease inhibitor 1 Hepatocyte growth factor activator inhibitor type 1 (HAI-1) SPINT1 HAI1
Observed MW:	58-70 kDa
Uniprot ID :	O43278

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded rat colon tissue slide using FNab10792(SPINT1 Antibody) at dilution of 1:200



HeLa cell lysates were subjected to SDS PAGE followed by western blot with FNab10792(SPINT1 Antibody) at dilution of 1:1000