

## TMEM41B antibody

### Product Information


Catalog No.:	FNab10806
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	transmembrane protein 41B
Synonyms:	Transmembrane protein 41B Protein stasimon TMEM41B KIAA0033
Observed MW:	22-28 kDa
Uniprot ID :	Q5BJD5

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:1000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse testis using FNab10806(TMEM41B antibody) at dilution of 1:100



Hela cell lysates were subjected to SDS PAGE followed by western blot with FNab10806(TMEM41B antibody) at dilution of 1:600