

progesterone antibody

Product Information

Catalog No.:	FNab10809
Size:	100µg
Form:	liquid
Purification:	protein A purified
Purity:	≥95% as determined by SDS-PAGE
Host:	mouse
Clonality:	monoclonal
Clone ID:	1G8
IsoType:	IgG1
Storage:	PBS pH 7.4, 2 - 8°C for short time, -20°C for 12 months.

Background

None

Immunogen information

Immunogen:	Progesterone-11-BSA
Synonyms:	Progesterone prog PRG PGN
Observed MW:	None
Uniprot ID :	None

Application

Reactivity:	universal
Tested Application:	ELISA, Lateral Flow, CLIA