

## TSG101 antibody

### Product Information


Catalog No.:	FNab10812
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	mouse
Clonality:	monoclonal
Clone ID:	1C6
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	tumor susceptibility gene 101
Synonyms:	Tumor susceptibility gene 101 protein ESCRT-I complex subunit TSG101 TSG101
Observed MW:	44 kDa
Uniprot ID :	Q99816

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC,FC
Recommended dilution:	WB: 1:1000 - 1:10000; IHC: 1:200 - 1:1000
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab10812(TSG101 antibody) at dilution of 1:500



Various lysates were subjected to SDS PAGE followed by western blot with FNab10812(TSG101 antibody) at dilution of 1:5000