

## LH antibody

### Product Information

Catalog No.:	FNab10817
Size:	1mg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3H3
IsoType:	IgG1
Storage:	PBS with 0.1% proclin300, pH 7.4, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	purified luteinizing hormone
Synonyms:	Lutropin subunit beta (Lutropin beta chain) Luteinizing hormone subunit beta (LH-B, LSH-B, LSH-beta) LHB
Observed MW:	None
Uniprot ID :	P01229

### Application

Reactivity:	Human
Tested Application:	ELISA, CLIA, GICA(detection antibody)
Recommended dilution:	Each laboratory should determine an optimum working titer for use in its particular application.