

TRIM72 antibody

Product Information


Catalog No.:	FNab10825
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	tripartite motif-containing 72
Synonyms:	Tripartite motif-containing protein 72 Mitsugumin-53 (Mg53) TRIM72 MG53
Observed MW:	53 kDa
Uniprot ID :	Q6ZMU5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:1000-1:5000; IP: 1:200-1:1000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsil tissue using FNab10825(TRIM72 antibody) at dilution of 1:100



Rat heart tissue were subjected to SDS PAGE followed by western blot with FNab10825(TRIM72 antibody) at dilution of 1:2000