

PDGFRA antibody

Product Information

Catalog No.:	FNab10830
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Platelet-derived growth factor receptor alpha
Synonyms:	Platelet-derived growth factor receptor alpha (PDGF-R-alpha, PDGFR-alpha) Alpha platelet-derived growth factor receptor Alpha-type platelet-derived growth factor receptor CD140 antigen-like family member A CD140a antigen Platelet-derived growth factor alpha receptor Platelet-derived growth factor receptor 2 (PDGFR-2) PDGFRA PDGFR2 RHEPDGFRA
Observed MW:	150 kDa
Uniprot ID :	P16234

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB : 1:1000-1:10000
Image:	



U251 cells were subjected to SDS PAGE followed by western blot with FNab10830(PDGFRα antibody) at dilution of 1:2000