

25-OH VD antibody

Product Information

Catalog No.:	FNab10834
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	sheep
Clonality:	monoclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.1% Sodium azide at 7.2-7.4, 2-8°C for 12 months.

Immunogen information

Immunogen:	25-OH VD
------------	----------

Application

Reactivity:	25-hydroxy vitamin D2/D3
Tested Application:	ELISA