

## PDX1 antibody

### Product Information


Catalog No.:	FNab10837
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Pancreas/duodenum homeobox protein 1
Synonyms:	Pancreas/duodenum homeobox protein 1 (PDX-1) Glucose-sensitive factor (GSF) Insulin promoter factor 1 (IPF-1) Insulin upstream factor 1 (IUF-1) Islet/duodenum homeobox-1 (IDX-1) Somatostatin-transactivating factor 1 (STF-1) PDX1 IPF1 STF1
Observed MW:	40 kDa
Uniprot ID :	P52945

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10837(PDX1 antibody) at dilution of 1:1000