

## ICA1 antibody

### Product Information


Catalog No.:	FNab10838
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Islet cell autoantigen 1
Synonyms:	Islet cell autoantigen 1 69 kDa islet cell autoantigen (ICA69) Islet cell autoantigen p69 (ICAp69, p69) ICA1
Observed MW:	70 kDa
Uniprot ID :	Q05084

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



U-87MG cells were subjected to SDS PAGE followed by western blot with FNab10838(ICA1 antibody) at dilution of 1:1000