

## FGF9 antibody

### Product Information


Catalog No.:	FNab10839
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Fibroblast growth factor 9
Synonyms:	Fibroblast growth factor 9 (FGF-9) Glia-activating factor (GAF) Heparin-binding growth factor 9 (HBGF-9) FGF9
Observed MW:	23 kDa
Uniprot ID :	P31371

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500 - 1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10839(FGF9 antibody) at dilution of 1:1000