

## NOTCH2 antibody

### Product Information


Catalog No.:	FNab10843
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	Notch homolog 2
Synonyms:	Neurogenic locus notch homolog protein 2 (Notch 2, hN2) Notch 2 extracellular truncation (N2ECD) Notch 2 intracellular domain (N2ICD) NOTCH2
Observed MW:	110 kDa
Uniprot ID :	Q04721

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab10843(NOTCH2 Antibody) at dilution of 1:50



HeLa cell lysates were subjected to SDS PAGE followed by western blot with FNab10843(NOTCH2 antibody) at dilution of 1:1000