

PYCARD antibody

Product Information

Catalog No.:	FNab10850
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9C1
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	PYD and CARD domain containing
Synonyms:	Apoptosis-associated speck-like protein containing a CARD (hASC) Caspase recruitment domain-containing protein 5 PYD and CARD domain-containing protein Target of methylation-induced silencing 1 PYCARD ASC CARD5 TMS1
Observed MW:	25 kDa
Uniprot ID :	Q9ULZ3

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab10850(ASC antibody) at dilution of 1:1000