

PMS1 antibody

Product Information

Catalog No.:	FNab10857
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probably involved in the repair of mismatches in DNA.

Immunogen information

Immunogen:	PMS1 postmeiotic segregation increased 1
Synonyms:	PMS1 protein homolog 1 DNA mismatch repair protein PMS1 PMS1 PMSL1
Observed MW:	60-100 kDa
Uniprot ID :	P54277

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500 - 1:2000
Image:	



HL-60 cell lysates were subjected to SDS PAGE followed by western blot with FNab10857(PSM1 antibody) at dilution of 1:1000