

PXDN antibody

Product Information


Catalog No.:	FNab10858
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Peroxidasin homolog
Synonyms:	Peroxidasin homolog Melanoma-associated antigen MG50 Peroxidasin 1 (hsPxd01) Vascular peroxidase 1 p53-responsive gene 2 protein PXDN active fragment PXDN KIAA0230 MG50 PRG2 PXD01 VPO VPO1
Observed MW:	165 kDa
Uniprot ID :	Q92626

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat lung using FNab10858(PXDN antibody) at dilution of 1:50



Various cell lysates were subjected to SDS PAGE followed by western blot with FNab10858(PXDN antibody) at dilution of 1:1000