

FGF3 antibody

Product Information


Catalog No.:	FNab10863
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Fibroblast growth factor 3
Synonyms:	Fibroblast growth factor 3 (FGF-3) Heparin-binding growth factor 3 (HBGF-3) Proto-oncogene Int-2 FGF3 INT2
Observed MW:	27 kDa
Uniprot ID :	P11487

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



Various cell lysates were subjected to SDS PAGE followed by western blot with FNab10863(FGF3 antibody) at dilution of 1:1000