

TP73 antibody

Product Information


Catalog No.:	FNab10864
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Tumor protein p73
Synonyms:	Tumor protein p73 p53-like transcription factor p53-related protein TP73 P73
Observed MW:	48-80 kDa
Uniprot ID :	O15350

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human esophagus using FNab10864(TP73 antibody) at dilution of 1:50



HeLa cell lysates were subjected to SDS PAGE followed by western blot with FNab10864(TP73 antibody) at dilution of 1:1000