

PBX2 antibody

Product Information


Catalog No.:	FNab10869
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	pre-B-cell leukemia homeobox 2
Synonyms:	Pre-B-cell leukemia transcription factor 2 Homeobox protein PBX2 Protein G17 PBX2 G17
Observed MW:	45, 60 kDa
Uniprot ID :	P40425

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Various cell lysates were subjected to SDS PAGE followed by western blot with FNab10869(PBX2 antibody) at dilution of 1:1000