

GNA15 antibody

Product Information


Catalog No.:	FNab10871
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Guanine nucleotide-binding protein subunit alpha-15
Synonyms:	Guanine nucleotide-binding protein subunit alpha-15 (G alpha-15, G-protein subunit alpha-15) Epididymis tissue protein Li 17E Guanine nucleotide-binding protein subunit alpha-16 (G alpha-16, G-protein subunit alpha-16) GNA15 GNA16
Observed MW:	43 kDa
Uniprot ID :	P30679

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Hela cell lysates were subjected to SDS PAGE followed by western blot with FNab10871(GNA15 antibody) at dilution of 1:1000