

ABCB5 antibody

Product Information


Catalog No.:	FNab10874
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	ATP-binding cassette sub-family B member 5
Synonyms:	ATP-binding cassette sub-family B member 5 ABCB5 P-gp P-glycoprotein ABCB5 ABCB5
Observed MW:	None
Uniprot ID :	Q2M3G0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human melanoma tissue using FNab10874(ABCB5 antibody) at dilution of 1:100