

## VSIG4 antibody

### Product Information


Catalog No.:	FNab10876
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	V-set and immunoglobulin domain-containing protein 4
Synonyms:	V-set and immunoglobulin domain-containing protein 4 Protein Z39Ig VSIG4 CRIG Z39IG
Observed MW:	34 kDa
Uniprot ID :	Q9Y279

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab10876(VSIG4 Antibody) at dilution of 1:500