

PELI3 antibody

Product Information


Catalog No.:	FNab10877
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	E3 ubiquitin-protein ligase pellino homolog 3
Synonyms:	E3 ubiquitin-protein ligase pellino homolog 3 (Pellino-3) PELI3
Observed MW:	50 kDa
Uniprot ID :	Q8N2H9

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500 - 1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab10877(PELI3 Antibody) at dilution of 1:800