

## OPRM1 antibody

### Product Information


Catalog No.:	FNab10889
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	Mu Opioid Receptor
Synonyms:	Mu-type opioid receptor (M-OR-1, MOR-1) Mu opiate receptor Mu opioid receptor (MOP, hMOP) OPRM1 MOR1
Observed MW:	60-70 kDa
Uniprot ID :	P35372

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain tissue using FNab10888(MOR antibody) at dilution of 1:100



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab10889(MOR Antibody) at dilution of 1:1000