

SOX4 antibody

Product Information


Catalog No.:	FNab10892
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Transcription factor SOX-4
Synonyms:	Transcription factor SOX-4 SOX4
Observed MW:	47 kDa
Uniprot ID :	Q06945

Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



SKOV3 cells were subjected to SDS PAGE followed by western blot with FNab10892(SOX4 Antibody) at dilution of 1:1000