

GYPA antibody

Product Information

Catalog No.:	FNab10895
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	None
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	glycophorin A(MNS blood group)
Synonyms:	Glycophorin-A MN sialoglycoprotein PAS-2 Sialoglycoprotein alpha GYPA GPA
Observed MW:	40 kDa
Uniprot ID :	P02724

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:500-1:2000
Image:	



Immunohistochemistry of paraffin-embedded human spleen tissue using FNab10895(GYPA antibody) at dilution of 1:1000



Human blood were subjected to SDS PAGE followed by western blot with FNab10895(GYPA antibody) at dilution of 1:1000