

SLC31A1 antibody

Product Information


Catalog No.:	FNab10897
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	High affinity copper uptake protein 1
Synonyms:	High affinity copper uptake protein 1 Copper transporter 1 (hCTR1) Solute carrier family 31 member 1 Truncated CTR1 form SLC31A1 COPT1 CTR1
Observed MW:	30 kDa
Uniprot ID :	O15431

Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution:WB:	1:500-1:1000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10897(SLC31A1 Antibody) at dilution of 1:600