

SOAT1 antibody

Product Information


Catalog No.:	FNab10902
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Sterol O-acyltransferase 1
Synonyms:	Sterol O-acyltransferase 1 Acyl-coenzyme A:cholesterol acyltransferase 1 (ACAT-1) Cholesterol acyltransferase 1 SOAT1 ACACT ACACT1 ACAT ACAT1 SOAT STAT
Observed MW:	55 kDa
Uniprot ID :	P35610

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10902(SOAT1 Antibody) at dilution of 1:1000