

TIGIT antibody

Product Information


Catalog No.:	FNab10904
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.06% PC300 and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	T-cell immunoreceptor with Ig and ITIM domains
Synonyms:	T-cell immunoreceptor with Ig and ITIM domains V-set and immunoglobulin domain-containing protein 9 V-set and transmembrane domain-containing protein 3 TIGIT VSIG9 VSTM3
Observed MW:	30 kDa
Uniprot ID :	Q495A1

Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Human plasma were subjected to SDS PAGE followed by western blot with FNab10904(TIGIT Antibody) at dilution of 1:1000