

Phospho-NFkB1 (S337) antibody

Product Information


Catalog No.:	FNab10905
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.06% PC300 and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	A phospho specific peptide corresponding to residues surrounding S337 of human NFkB1
Synonyms:	Nuclear factor NF-kappa-B p105 subunit DNA-binding factor KBF1 EBP-1 Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1 Nuclear factor NF-kappa-B p50 subunit NFkB1
Observed MW:	120 kDa
Uniprot ID :	P19838

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab10905 at dilution of 1:1000