

## HDAC7 antibody

### Product Information


Catalog No.:	FNab10906
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.06% PC300 and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	histone deacetylase 7
Synonyms:	Histone deacetylase 7 (HD7) Histone deacetylase 7A (HD7a) HDAC7 HDAC7A
Observed MW:	120 kDa
Uniprot ID :	Q8WUI4

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10906(HDAC7 Antibody) at dilution of 1:1000