

## SLC26A3 antibody

### Product Information


Catalog No.:	FNab10912
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Solute carrier family 26 member 3
Synonyms:	Chloride anion exchanger Down-regulated in adenoma (Protein DRA) Solute carrier family 26 member 3 SLC26A3 DRA
Observed MW:	85 kDa
Uniprot ID :	P40879

### Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



NCI-H125 cells were subjected to SDS PAGE followed by western blot with FNab10913 at dilution of 1:1000.