

Phospho-RB (S780) antibody

Product Information


Catalog No.:	FNab10919
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	A synthetic phosphorylated peptide around S780 of human RB
Synonyms:	Retinoblastoma-associated protein p105-Rb p110-RB1 pRb (Rb) pp110 RB1
Observed MW:	110 kDa
Uniprot ID :	P06400

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab10919 at dilution of 1:1000. NIH/3T3 cells were treated by UV at room temperature for 15-30 minutes.