

LCN2 antibody

Product Information


Catalog No.:	FNab10922
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Neutrophil gelatinase-associated lipocalin
Synonyms:	Neutrophil gelatinase-associated lipocalin (NGAL) 25 kDa alpha-2-microglobulin-related subunit of MMP-9 Lipocalin-2 Oncogene 24p3 Siderocalin p25 LCN2 HNL NGAL
Observed MW:	22 kDa
Uniprot ID :	P80188

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10922(NGAL Antibody) at dilution of 1:1000