

## GGT1 antibody

### Product Information


Catalog No.:	FNab10923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Glutathione hydrolase 1 proenzyme
Synonyms:	Glutathione hydrolase 1 proenzyme Gamma-glutamyltransferase 1 Gamma-glutamyltranspeptidase 1 (GGT 1) Leukotriene-C4 hydrolase Glutathione hydrolase 1 heavy chain Glutathione hydrolase 1 light chain GGT1 GGT
Observed MW:	24 kDa, 40-70 kDa
Uniprot ID :	P19440

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10923(GTT1 Antibody) at dilution of 1:1000