

## anti- INSM1 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab10934   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | Insulinoma-associated protein 1             |
| Synonyms:    | INSM1, IA1, Insulinoma-associated protein 1 |
| Observed MW: | 53 kDa                                      |
| Uniprot ID : | Q01101                                      |

### Application

|                           |              |
|---------------------------|--------------|
| Reactivity:               | Human, Mouse |
| Tested Application:       | WB           |
| Recommended dilution: WB: | 1:500-1:2000 |
| Image:                    |              |


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Various lysates were subjected to SDS PAGE followed by western blot with FNab10934(INSM1 antibody) at dilution of 1:1000