

SLC1A5 antibody

Product Information


Catalog No.:	FNab10940
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Solute carrier family 1 member 5
Synonyms:	Neutral amino acid transporter B(0) (ATB(0)) Baboon M7 virus receptor RD114/simian type D retrovirus receptor Sodium-dependent neutral amino acid transporter type 2 Solute carrier family 1 member 5 SLC1A5 ASCT2 M7V1 RDR RDRC
Observed MW:	70 kDa
Uniprot ID :	Q15758

Application

Reactivity:	Human, Mouse
Tested Application:	WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab10940(ASCT2 Antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10940(ASCT2 Antibody) at dilution of 1:1000.