

CD200R1 antibody

Product Information

Catalog No.:	FNab10941
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Cell surface glycoprotein CD200 receptor 1
Synonyms:	Cell surface glycoprotein CD200 receptor 1 CD200 cell surface glycoprotein receptor Cell surface glycoprotein OX2 receptor 1 CD200R1 CD200R CRTR2 MOX2R OX2R
Observed MW:	45 kDa
Uniprot ID :	Q8TD46

Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:	WB: 1:500-1:2000
Image:	

Various lysates were subjected to SDS PAGE followed by western blot with FNab10941(CD200R Antibody) at dilution of 1:1000.

